



Guidelines for SAMSEG QC

- It is not necessary with a detailed QC of all the segmented regions, but check that the overall segmentation looks reasonable.
- In cases where SAMSEG has failed horribly, the reason is usually misregistration of the image to the template. This can be caused by excessively large field-of-views or other issues with the MR image.
- Minor errors where non-brain is classified as brain or similar is not grounds for exclusion for our purposes as SAMSEG is only used for estimating sbTIV.
- The next slides show examples of acceptable and bad segmentations from a variety of scanners.
- Do not hesitate to contact Claudia or Stener if you have questions.


Acceptable



Acceptable



Bad



Bad

